2. Apply **E.D.T.A. Trissodium GEL** in dental wall and wait the product start acting during 2 to 3 minutes;
3. Repeat procedure, if necessary;
4. Rinse thoroughly with physiologic solution;
5. Suture if necessary;
6. Rinse with physiologic solution or deionized water;
7. Apply surgical cement if necessary.

III - Prosthesis:

1. After canal prepare with adequate burs, place **E.D.T.A. Trissodium GEL** and leave it acting during 5 minutes;
2. Rinse with an irrigator solution of your choice (a.e. Milton's Liquid, Dakin's Liquid, Chlorinated Soda, Tergental and Tergipol);
3. Repeat the procedure if necessary;
4. Cement with the product of your choice.

SPECIAL CARE:

Dispose the product in accordance with local regulations, ensuring its complete distortion, preventing its reuse and environment damages.

CONTRAINDICATION:

E.D.T.A. Trissodium GEL should not be used in patients is know to have sensibility to any of the ingredients listed.


WARNINGS:

Do not use material after expiry date.

The manufacturer is not liable for damages caused by improper use or use not provided at these instructions.

"KEEP MATERIAL OUT OF THE REACH OF CHILDREN".

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COMPOSITION:

Ethylene diamine tetraacetic Acid; Sodium Hydroxide; Deionized Water and Thickener.

INDICATIONS:

To facilitate endodontic instrumentation;
In the etching of the tooth root in periodontic surgery.
As an auxiliary in the root scaling and to help the cleaning of the conduit for the cementing of pins.

TECHNICAL INFORMATION:

The Trissodium E.D.T.A. is a widely used product in the preparation of root canal walls, previously to the filling. It acts demineralizing dental particles by chelating the ions of calcium and magnesium, facilitating its dissolution and absorption, revealing itself as an effective chelating and lubricant agent. It is also effective in the dental scrap removal during the periodontal therapy, displaying the collagen and facilitating the tack of the conjunctive tissue treated on the root surface. This composition decalcifies the dentine to a depth of 10 to 30 µm in 5 minutes. Its exclusive formulation in gel form and its high concentration (24%) results in a bigger effectiveness of the product. The 24% concentration does not promote inconvenience for the dentine or the periapex, but it confers to the product a faster and efficient action. The product takes from 2 to 3 minutes to start making effect and for 5 to 6 minutes it reacts with the dentine and is saturated with calcium ions, losing its effect and neutralizes itself. The pH of the product is approximately 8,5 to 9,0. The viscosity of the gel was especially developed to provide the draining of the product for all the wall of the canal root. The gel is hydro-soluble allowing its complete elimination of the canal root through the cleanness with the indicated irrigator solutions.

USE INSTRUCTIONS:

I- Endodontics:

1. Place a small quantity of **E.D.T.A. Trissodium GEL** in the entrance of radicular conduit;
2. Wait until the product start acting for 2 to 3 minutes and start instrumentation. After 5 to 6 minutes the product is going to lose it's effect, because it already reacted with the dentin and was neutralized;
3. Rinse and if necessary repeat the procedure as many time as necessary, placing another quantity of the product;
4. Keep root canal moistened with **E.D.T.A. Trissodium GEL** during all the operating procedure;
5. After finishing instrumentation, rinse thoroughly with an irrigator solution of your choice. (Sodium hypochlorite solution, Tergipol or Tergental).

II - Periodontics:

1. Remove granulation tissue;