INDICADOR BIOLÓGICO 8 E 24 HORAS

INDICADOR BIOLÓGICO 8 E 24 HORAS / 8 AND 24-HOUR BIOLOGICAL INDICATOR

Registro de Controles Biológicos

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REFERÊNCIAS:

ISO 11138-1:2017 - Sterilization of health care products - Biological indicators - Part 1: General requirements.
ISO 11138-3:2017 - Sterilization of health care products - Biological indicators - Part 3: Biological indicators for moist heat sterilization processes.

8 hours and 24 hours Biological Indicator

EN

PRESENTATION:

Box with 10 plastic tubes with paper strips containing *Geobacillus stearothermophilus ATCC® 7953™* spores and a sealed glass ampoule with culture medium.

INDICATION:

Recommended for periodic monitoring of sterilization cycles and validation of steam autoclaves.

DESCRIPTION:

The ALLPRIME BIOLOGICAL INDICATOR serves the purpose of monitoring steam sterilization in healthcare facilities. It contains heat-resistant bacterial spores (*Geobacillus stearothermophilus ATCC® 7953™*), offering greater reliability in the sterilization process.

COMPOSITION:

· Purple culture medium, contained in the glass vial · *Geobacillus stearothermophilus ATCC® 7953™* with a minimum population of 105 to 106 spores inoculated on paper strips· Adhesive label with a chemical indicator (thermochromatic strip).

DIRECTIONS FOR USE:

1. Take one Biological Indicator tube from the box for each cycle or autoclave you wish to test (test tube). Also, take one control tube for each incubation cycle you perform and/or incubator you use. Identify the equipment, cycle, and date on each label. · Test tube: Place the tube, with the intact ampoule, inside a sterilization pouch made of surgical-grade paper for subsequent autoclave sterilization. · Control tube: For comparative purposes, keep the control tube outside the autoclave, with the ampoule also intact, and at room temperature. 2. Insert the sterilization pouch with the test tube into the most critical location in the autoclave (usually, the upper and/or near the door) along with the materials to be sterilized, and follow the manufacturer's instructions for starting the sterilization cycle. 3. After the cycle is complete, remove the sterilization pouch from the autoclave and let it cool for at least 15 minutes before proceeding. 4. Ensure that there has been a color change in the chemical indicator on the label, from pink to brown, after the sterilization process. 5. Remove the test tube from the sterilization pouch and activate the ampoules (test and control) as follows: · Apply enough pressure to the tube to break only the internal glass ampoule and release the culture medium. When using a mini-incubator, you can achieve the breakage by compressing the tubes against the internal walls of the incubator holes. Alternatively, you can use an ampoule breaker. · Be careful not to break the plastic tube and always keep the tube in an upright position. Ensure you bring the culture medium into contact with the paper strips containing the spores. If necessary, tap lightly with your finger to remove air bubbles and ensure proper contact. 6. After breaking, start the incubation process of the ampoules for the time indicated on the label, at a controlled temperature of 55 °C to 60 °C. 7. Highlight the labels and record the results on the control sheet, located on the back of this manual, following the demonstrated model.

INTERPRETATION OF RESULTS (BIOLOGICAL INDICATOR FOR 8H AND 24H):

Test Tube – The ampoule that has undergone autoclave sterilization should remain purple after incubation, and the label strip should have a brown hue.

· The persistence of the purple color in the test sample indicates the absence of bacterial growth, which is expected. In this case, spore destruction has occurred, confirming adequate sterilization, according to Results Table (a) and (b).

· The change in color to a yellow hue, which is not expected in this situation, as it suggests inadequate sterilization, according to Results Table (c) and (d). Action: In this case, you should stop using the autoclave, analyze the physical parameters and possible equipment oversignals bacterial growths, and repeat the sterilization cycle with a new ampoule. If the yellow hue persists, we recommend sending the equipment for technical assistance. After repairing it, you should repeat the sterilization procedure.

Control Tube – The ampoule that has not undergone the autoclave sterilization process should change color to a yellowish hue after incubation, and the label strip should remain pink.

· A color change to a yellowish hue signals bacterial growth, which usually occurs before the incubation time described on the label. This is the expected result for the control tube, according to Results Table (a) and (c).

· If, after incubation of the control tube, there is no color change to yellow, the result will be negative for bacterial growth, invalidating the test, as sterilization cannot be confirmed, according to Results Table (b) and (d). Action: This lack of color change may be due to a malfunction of the incubator, improper breakage of the glass ampoule, or incorrect product storage. Also, make sure the control tube has not been sterilized. It is important to investigate the cause of this result and take corrective measures.

RESULT:

The test will be considered approved if the test vial (sterilized) shows a purple hue and the control vial shows a yellow hue. Only situation (a) in the Results Table represents a satisfactory result for both viais. Any result that does not meet these criteria will be considered failed.

INTERPRETATION OF RESULTS (SPECIFICS FOR 8-HOUR BIOLOGICAL INDICATOR):

1. The quick detection of biological indicators, with readings within 8 hours (ALLPRIME 8h Biological Indicator), occurs through the gradual change in color of the growth medium, as represented in the Color Reference Table presented below, in this manual. 2. The color verification should be performed under white light and can be initiated from the second hour of incubation. 3. Confirmation of the negative result for bacterial growth should occur from the 7th hour of incubation, noting the continued presence of the purple color (expected for the test tube). 4. The positive result for bacterial growth is indicated by the change in the medium to yellow and should be observed within a maximum of 4 hours after the start of incubation, remaining in this color until completing 8 hours (expected for the control tube). 5. It is important to emphasize that the tubes should not be removed from the incubator before completing the total incubation time to avoid interference in the results.

PRECAUTIONS:

· It is essential to observe and strictly follow the provided instructions for use. · Access the product page on the website <https://allprimedental.com> to access illustrated and video instructions. · It is advisable to use only one control tube per incubation cycle, even if multiple test tubes are used. · Ensure that the control ampoules are from the same batch as the test ampoules. · Turn on the incubator at least one hour before use to ensure it has reached a temperature of 55 °C a 60 °C by the time of incubation. · The incubation process must be initiated within a maximum of two hours after autoclave sterilization. · On the ampoule label, check for a change in the chemical indicator from pink to brown after the sterilization process. The color change on the label only confirms exposure to steam at 121°C or above and does not indicate the conformity of the sterilization process. · If the liquid evaporates during the sterilization or incubation process, disregard the results and repeat the operation. · In case of recurrent failed results, contact the biological indicator manufacturer.

WARNINGS:

· The ampoules should not be reused and should be disposed of properly after autoclave sterilization. · Keep the material out of reach of children. · This material has been manufactured for the specific intended use and must be handled according to the instructions for use. · The manufacturer is not responsible for damages caused by other uses or incorrect handling not described in this manual. · Additionally, the user is obligated to verify, before use and under their responsibility, whether this material is compatible with the intended use, especially when such use is not indicated in the provided instructions. · Do not expose the material directly to sunlight. · Do not use the material if the packaging is damaged. · This product should be used only in accordance with the instructions for use. · Use personal protective equipment (PPE) when handling the product, including gloves and protective eyewear.

STORAGE AND TRANSPORT CONDITIONS:

Store the product properly, avoiding direct exposure to light, maintaining it within a temperature range of 15 °C to 30 °C, and with a relative humidity between 35% and 60%.

DISPOSAL METHOD:

You must sterilize all materials used in microbiological control for a minimum of 30 minutes at a temperature of 121 °C and then properly discard them.

BIOLOGICAL INDICATOR RESISTANCE DATA FOR STEAM STERILIZATION:

Population (average number of spores/unit), survival time in minutes, death time in minutes, and D-value (121 °C) are determined for each batch. This information appears on the adhesive label, present in the product box

REFERENCES:

ISO 11138-1:2017 - Sterilization of health care products - Biological indicators - Part 1: General requirements.
ISO 11138-3:2017 - Sterilization of health care products - Biological indicators - Part 3: Biological indicators for moist heat sterilization processes.

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