Perfil

PORTUGUÊS

1. Descrição

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ENGLISH

1. Product Description

Perfil is a product based on condensation silicone indicated for precision dental impression or edentulous arch, aiming to elaborate dental prothesis models total or partial. Perfil is presented in the versions: Denso (High Consistency), Light Body (Low Consistency) and Cub (Universal Catalyst).

1.1. Indications

- Double step technique
- Single step technique
- Technique with copper ring
- Impressions for study models, orthodontic models, dies, protective isolation during staining of the teeth of the prosthesis in the oven.

1.2. Composition

- Perfil Denso: Polydimethylsiloxane, Mineral Oil, Inorganic Filler, Pigment and Flavor
- Perfil Fluido: Polydimethylsiloxane, Inorganic Filler, Pigment and Flavor
- Perfil Cub: Hardener for Condensation Silicone, Bis((2-ethyl-1-oxohexyl)oxy) diocetylstannane, Vaseline, Silicon Dioxide and Pigment.

1.3. Presentation

Perfil Kit:

- 01 jar of Perfil Denso with 1 Kg
- 01 tube of Perfil Fluido with 120 g
- 01 tube of Perfil Cub with 50 g
- 01 measuring spoon
- 01 Instructions for Use

Perfil Fluido Refill:

- 01 tube of Perfil Fluido with 120 g
- 01 Instructions for Use

Perfil Denso Refill:

- 01 jar of Perfil Denso with 1 Kg
- 01 measuring spoon
- 01 Instructions for Use

Perfil Cub Refill:

- 01 tube of Perfil Cub with 50 g

01 Instructions for Use

1.4. Flavor: mint

2.1. Instructions

2.1.1. Double Impression

- Select a rigid tray, with the capacity of providing 3-4 mm of thickness of impression material.
- Clean and dry the field, if required moving the gum away using retraction cords.
- Using a measuring spoon, remove the desired quantity of Perfil Denso from the jar, ensuring that the measuring spoon is properly filled without empty spaces or excess of material.

- Press the Perfil Denso against the palm of your hand and mark the diameter of the measuring spoon against the Perfil Denso.
- For each spoonful of Denso, dispense the catalyst (Perfil Cub) in the proportion equivalent to the diameter marked.
- Mix the Perfil Denso and Perfil Cub with your fingertips to obtain a uniform coloring. The time required to execute the complete homogenization is 30 seconds.

- Accommodate the mixture into the tray and position the tray on the area to be molded, exerting a pressure for 4-5 seconds and maintaining it stable for another 3 minutes.
- Remove the tray from the oral cavity.

NOTE: Before making the silicone molding dense is necessary to promote a relief with PVC or other plastic. The relief may be performed after molding using a cutting instrument.

- Dispense on a glass plate the required quantity of Perfil Fluido, being careful that the material is dispensed in a straight line.
- For each mark of Perfil Fluido, dispense the same quantity of Perfil Cub. The use of a larger quantity of catalyst accelerates the setting time.

- Incorporate the fluid paste into the catalyst mixing vigorously with a spatula so as to homogenize all the material in 30 seconds.

- Using an impression syringe, take part of the material to the cavity preparation. The rest of the material must be dispensed inside the preliminary Denso impression.

- Position the tray in the place desired exerting a pressure for 4-5 seconds and maintains it stable for another 4 minutes before removing it from the oral cavity.

- After removing the impression, wash it with water and disinfect it before pouring the model.

2.2. Simultaneous Molding

- Select a rigid tray, with the capacity of pro-

viding 3-4 mm of thickness of impression material.

- Clean and dry the field, if required moving the gum away using retraction cords.
- Using a measuring spoon, remove the desired quantity of Perfil Denso from the jar, ensuring that the measuring spoon is properly filled without empty spaces or excess of material.
- Press the Perfil Denso against the palm of your hand and mark the diameter of the measuring spoon against the Perfil Denso.
- For each spoonful of Denso, dispense the catalyst (Perfil Cub) in the proportion equivalent to the diameter marked.
- Dispense on a glass plate the required quantity of Perfil Fluido, being careful that the material is dispensed in a straight line.
- For each mark of Perfil Fluido, dispense the same quantity of Perfil Cub. The use of a larger quantity of catalyst accelerates the setting time.

- Mix the Perfil Denso and Perfil Cub with your fingertips to obtain a uniform coloring. The time required to execute the complete homogenization is 30 seconds.
- Accommodate the mixture into the tray.
- Simultaneously, the assistant or professional must start the mixing of Perfil Fluido with com o Perfil Cub stirring vigorously with a spatula so as to homogenize all the material in at most 30 seconds.

- Using an impression syringe, inject the mixture around the cavity preparation.
- Position the tray in the place desired exerting a pressure for 4 thru 5 seconds and maintains it stable for another 4 minutes before removing it from the oral cavity.
- After removing the impression, wash it with water and disinfect it before pouring the model.

Type	Perfil Denso
	Putty Consistency
Mixture time	0'30"
Working time	1'30"
Time in the mouth	3'00"
Setting time	4'30"

Type	Perfil Fluido
	Light Consistency
Mixture time	0'30"
Working time	1'30"
Time in the mouth	4'00"
Setting time	5'30"

2.3. Desinfection

After the impression had being rinsed with water, any dental disinfectant chemical solution suitable commercial can be used.

3. Precautions, Restrictions and Warnings

3.1. Precautions

- All dental material should only be used by dentists and authorized dental labs.
- The setting time is considered for a temperature 23°C. An increased temperature accelerates the setting time and a decreased temperature reduces the setting time.
- An increased quantity of Perfil Cub reduces the working time.
- Avoid the contact eyes. It is advisable to use safety goggles for patients, assistants and professionals during the all clinical procedure.
- Perfil Cub is irritating.
- Keep the product away from children.
- Single use.

3.2. Restrictions

There are no contraindications when the product is used by a certified professional following the Instructions for Use. Avoid using the product on any patients with an allergic history to one or more components of the product.

4. Expiration Date

Do not use after the expiration date. Under storage conditions, Perfil has a shelf life of 3 years.

5. Storage

Store at room temperature, in a dry place and protected from light and heat. Replace cap tightly after use. Inadequate storage conditions will shorten the shelf life and may lead to a malfunction of the product.

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