## MAXDAM - Barreira gengival fotopolimerizável

**DESCRIÇÃO:**

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**DESCRIPTION:**

Hydrogen peroxide or the whitening gel, when in contact with the gums, can damage the tissues and cause great discomfort, thus interfering with the whitening result. Maxdam acts to avoid this discomfort, guaranteeing a better result in the teeth whitening treatment. It is a gel-based product that acts as a light-curing gingival barrier and replaces the rubber sheet when gum isolation is relative. The main purpose of the product is to protect the tissues from contact with hydrogen peroxide or whitening gel. It is suitable for isolating the sulci and gums during in-office tooth whitening procedures.

**PRESENTATIONS:**

Maxdam Blue: 1 syringe with 2 g and 3 application tips.

Maxdam Green: 1 syringe with 2 g and 3 application tips.

**INDICATION:**

It is indicated for the protection of gum tissue, replacing the rubber sheet in cases where the clinical procedure may cause some discomfort or irritation due to contact with the products used.

**COMPOSITION:**

HEMA, silica, urethane dimethacrylate monomers, inert filler, pigments and photoinitiators.

**MODE OF USE:**

- Choose the color of Maxdam according to your preference.
- Check that the packaging has not been damaged, if it has, discard it.
- Aseptinize the area to be worked on, drying the area of the gums and teeth that will be protected.
- Use a lip retractor.
- Perform relative isolation, thoroughly drying the area of the gums and teeth that will be protected.
- Using the applicator tip, apply a layer approximately 1 to 1.5 mm thick and up to 3.5 mm wide of gingival barrier to the gums, including the previously dried papillae.
- The barrier should cover approximately 0.5 to 1 mm of the tooth surface for full coverage and patient protection.
- Light-cure the product immediately after application.
- Cure for 15 seconds for each group of 3 to 4 teeth. Scanning light curing can be used, as long as you use light cures with a power greater than 800 mW/cm<sup>2</sup>.
- Check that the entire gingival surface is protected and that the entire barrier has been light-cured properly.
- Perform whitening following the protocols of each manufacturer.
- To remove the gum barrier, simply push on it with a blunt instrument until it comes off easily in whole or in large pieces.

**WARNINGS:**

- For use by the dental surgeon only.
- Do not use the product after its expiration date.
- Do not reuse empty packaging.
- Keep out of reach of children.
- This material has been manufactured for dental use only and must be handled in accordance with the instructions for use. The manufacturer is not responsible for damage caused by other uses or by incorrect handling not described in this manual. In addition, the user is obliged to check, prior to use and

under their own responsibility, that this material is compatible with the intended use, especially when this use is not indicated in the present instructions for use.

- Maxdam should not be swallowed or aspirated. If swallowed or aspirated, seek medical attention immediately.
- In the case of scanning applications, use a curing light with a power greater than 800 mW/cm<sup>2</sup>. Photopolymerizers with lower power may not guarantee complete polymerization of the product.
- Strictly follow the procedures provided in the instructions for use.
- Do not use the product on the gums of patients who have recently undergone oral surgery or who have any kind of gum irritation.
- In the event of acute allergic reactions, discontinue use of the product.
- Discard the tip after use.
- Incomplete polymerization of the product or excessive time may lead to slight gum irritation. The dental surgeon should monitor the patient and pay attention to their sensitivity.
- Do not handle the Maxdam syringe with the light of the reflector directly focused on the open packaging.
- Do not project the light from the light curing device directly onto the open packaging.

**PRECAUTIONS:**

- Wear PPE when handling the product, including gloves and goggles. It is recommended that patients also wear goggles.
- Do not allow loose polymerized objects and parts to be swallowed or aspirated.
- Contact of unpolymerized Maxdam with skin, mucous membranes and eyes should be avoided. When not yet polymerized, the product may have a slightly irritating effect.
- If the product comes into contact with your hands, wash with plenty of water.
- In the event of allergic reactions, discontinue use immediately. Seek medical advice if necessary.
- Product for single use only. After polymerization, the product does not return to its original consistency.
- No reprocessing.
- Do not expose the material to direct sunlight.
- Replace cap on syringe immediately after use.

**CONTRAINDICATIONS:**

Do not use this material in patients who are allergic or hypersensitive to any component of the formula.

**ADVERSE EFFECTS:**

During the clinical procedure, the patient may feel a slight variation in the temperature of the product. The sensation of heat and burning may increase the longer the photoactivation time. This effect is common with resin products based on methacrylic monomers, due to the polymerization and contraction of the resin. In such cases, the dentist should carry out scanning polymerization in order to reduce the heating effect. Maxdam is composed of methacrylic monomers, and some patients may be sensitive to this class of compounds. The dental surgeon must assess the sensitivity of each patient. In cases of incomplete polymerization, whitening gel may pass between the barrier and the gums, causing damage to the patient. The dental surgeon must carry out the light-curing process properly and ensure that the barrier is perfectly fixed before starting to apply the whitening agent.

**STORAGE CONDITIONS:**

Always store the product in its original, unopened packaging at room temperature. Protect from direct sunlight. Do not store near products containing eugenol.

**TRANSPORT CONDITIONS:**

The product can be transported between -28 °C and 55 °C for up to 15 days without compromising its initial performance.

**VALIDITY:**

If packaged properly, the product has a shelf life of 24 months from the date of manufacture. Do not use after the expiration date.

**DISPOSAL MODE:**

Dispose of in accordance with local legislation, ensuring its complete de-characterization in order to prevent its reuse and environmental damage.

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Maquira Indústria de Produtos Odontológicos S.A.