INSTRUCTIONS FOR USE
MIR - MICRODONT INTERPROXIMAL REDUCTION
MIR - MICRODONT REDUCCION INTERPROXIMAL
MIR - MICRODONT REDUÇÃO INTERPROXIMAL

MS: 80676920040

ENGLISH

Indication:

The MIR is indicated mainly for orthodontic stripings (technique in which it is desired to create interproximal space through enamel wear) and can also be used for interproximal finishing and for removal of excess cements and restorative materials. It has a sequence of six arcs, one saw and five abrasive strips ranging from extrafine to extra thick for MIR 1.0 and five arcs, being a saw and four abrasive strips ranging from extrafine to thick for MIR 2.0.

Fundamentals of its operation and action:

Abrasion

Instructions for use:

The MIR was developed to assist in the interproximal opening of dental elements for orthodontic stripping. Each saw arc or abrasive strip has a numbering (in millimeters) corresponding to the size of the opening, to obtain the said opening must be used following the ascending order of arcs:

MIR 1.0:

- White – Serrated – 0,05mm
- Yellow – Extra Thin – 0,10mm
- Red – Thin – 0,20mm
- Blue – Regular – 0,30mm
- Green – Thick – 0,40mm
- Lilac – Extra Thick – 0,50mm
- The MIR 1.0 arches can be used on their own handle, with the aid of the cable or in counter-angle.

MIR 2.0:

- White – Serrated – 0,05mm
- Yellow – Extra Thin – 0.05mm single and 0.10mm double face
- Red – Thin – 0.10mm single face and 0.20mm double face
- Blue – Regular – 0.15mm single face and 0.30mm double face
- Green – Thick – 0.20mm single face and 0.40mm double face
- MIR 2.0 arcs are exclusively manual.

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
Sequence of Use:

Manual cleaning procedure:

- Use a reserved area specifically to perform cleaning.
- Rinse the MIR in running water and keep them immersed in a bath of enzymatic detergent solution, following the manufacturer's instructions.
- After soaking and keeping it immersed, brush well away from the body with the neutral cleaning agent. Care should be taken to avoid spreading contaminants from the spray during brushing procedures. Do not use wire brushes.
- Special care should be taken to clean the cracks and other hard-to-reach areas. Inspect visually to confirm debris removal. Repeat the cycle if necessary.
- Rinse the MIR thoroughly under running water until it is visibly clean.
- Dry the device with clean compressed air or paper towel.

Ultrasonic cleaning procedure

- Prepare an enzymatic detergent solution by following the agent manufacturer's instructions for the correct concentration, exposure time, temperature and water quality.
- Place the MIR in the proper place of the ultrasonic equipment, following the instructions for use of the equipment manufacturer
- Completely immerse the device in the cleaning solution and program the equipment for at least fifteen (15) minutes.
- Rinse thoroughly under running water until it is visibly clean.
- Dry the device with clean compressed air or paper towel.

Sterilization requirements:

In case of reuse it is recommended to sterilize the product in pressure saturated steam autoclave at 121 ° C for 20 minutes, properly packed in its own autoclave packaging, or the sterilization parameters must be followed in accordance with the relevant policies / rules of the country where the device is being used. It is the sole responsibility of the professional to keep the autoclave equipment updated to ensure the correct sterilization of the product.

Contraindication:

- Cleaning agents with chlorine or chloride as the active ingredient are corrosive to stainless steel and should not be used.
- Do not use cold sterilization methods (chemical process) for the sterilization of the MIR. These chemical agents can cause corrosion of the products.

Warnings and Precautions:

- Before use, inspect if the product is in perfect condition of use;
- Non-sterile product. 3R recommends single use.
- Sterilize the product before use.
- During the operation the professional must use all personal protective equipment, including goggles.
- Wear goggles on the patient.
- The Product is considered to be a contaminated instrument and proper handling precautions are required. Wear PPE, including gloves, masks, safety goggles, and a long-sleeved lab coat when handling, just as is done with all other dental instruments.

INSTRUCTIONS FOR USE
MIR - MICRODONT INTERPROXIMAL REDUCTION
MIR - MICRODONT REDUCCION INTERPROXIMAL
MIR - MICRODONT REDUÇÃO INTERPROXIMAL

MS: 80676920040

- LOT Number: It is important to keep the batch number, identified on the product label, to be informed in any occurrences with the manufacturer.
- The disposal of the material must follow the relevant local legislation.
- Keep the package in a dry place, away from light and away from the weather. Keep the products in their original packaging until the moment of their use.

Conservation and storage:




- The product must be stored in an airy place, dry and protected from light, in order to maintain the ideal conditions of storage and transport, as well as its physical and chemical integrity.
- Keep the product in its original packaging until the time of use. Inspect the packaging before use.
- The operator must be informed of the content. The care to be taken in the transport is demonstrated in the packaging, thus ensuring the protection of the instrument from shipment to delivery to the customer. Packets should not be left to the action of time.

Other Additional Information:

- The sterilization processes described above have been validated and are adequate and effective for the use of MIR. 3R is not responsible if the instructions for use are not followed properly, generating unwanted results in the efficacy and safety of the product/patient.
- If necessary, 3R provides the instructions for use in printed format, at no additional cost, upon request by e-mail (administrativo@3r.ind.br) or by phone: (+55 11-5525-0590 / +55 0800-773-5525), and we will return within a maximum of 7 working days from the request.
- Always check on the site: www.microdont.com.br, the most up-to-date version of the usage instruction corresponding to the product.

Composition: Stainless steel, diamond and polycarbonate.

Symbology:

Symbol	Symbol Title	Description
	Expiration / Use-by date	Indicates the date after which the health product should not be used.
	Consult instructions for use	Indicates the user's need to refer to the instructions for use.
	Date of manufacture	Indicates the date when the health product was manufactured
	Manufacturer	Indicates the Manufacturer of the product for the health as defined in the Directive of the European Community 90/385EEC, 93/42/EEC and 98/79/EC.
	Batch Code	Indicates the manufacturer's batch code so that the batch or departure can be identified.
	Catalog number	Indicates the manufacturer's reference number so that the product for the health can be identified.
	Sterilized using irradiation	Indicates that the product is sterilized
	Non-Sterile	Indicates that the product is not sterile
	Caution	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself
	Authorized representative in the European Community	Indicates the Authorized representative in the European Community

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This review cancels and replaces the previous revision.

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