VOCO

Amaris® Gingiva

CE 0482

EN

Instructions for use

MD EU Medical device

In accordance with DIN EN ISO 4049

Product description:

Amaris Gingiva is a light-curing radiopaque restorative which contains 80% w/w fillers in a methacrylate matrix and is cured with blue light. **Amaris Gingiva** is available in syringes. **Amaris Gingiva** is used with a dentine / enamel bond.

Amaris Gingiva is used to treat cervical defects and exposed, discoloured or hypersensitive necks of teeth, especially in the visible anterior area. The gingival-like shade, together with the **Opaker** system, allows the correct reproduction of the gingiva also in its shade. The system comprises of three **Opakers** in different shades and a **gingiva-hued composite restorative in a basic shade**.

The light-curing **Opakers** are used to mask discoloured cervical areas and to individually adapt the restoration to the surrounding gingiva.

The light-curing **Amaris Gingiva** allows for a simplified and fast restorative technique.

Indications:

- Class V restorations (cervical caries, root erosion, V-shaped defects) especially in cases of retreating gingiva e.g. as a result of periodontal disease.
- Covering of exposed, discoloured or hypersensitive necks of teeth, especially in the visible anterior area.
- Highly aesthetic correction of malpositions of teeth, for direct veneers, facings and corrections of the red-white aesthetics.

Contraindications:

Amaris Gingiva contains methacrylate and BHT. **Amaris Gingiva** should therefore not be used in patients with a known hypersensitivity (allergy) to these ingredients.

Patient target group:

Amaris Gingiva is suitable for use in all patients without any age or gender restrictions.

Performance features:

The product's performance features satisfy the requirements of the intended use and the relevant product standards.

Use:

Amaris Gingiva should only be used by a professionally trained dental practitioner.

Use:

Preparation:

Clean the teeth being treated. Ensure that the working field is sufficiently dry. Use of a rubber dam is recommended.

Cavity preparation:

As a matter of principle, minimally invasive preparation (adhesive restorative technique) should be carried out to conserve healthy tooth substance. It is not necessary to prepare caries-free cervical defects. Here, cleaning with a cleaning paste or pumice and a rubber cup or rotating brushes is sufficient.

The use of a retraction cord can be helpful.

Cleaning:

- Remove all residues from the cavity with water spray.
- Any contamination with blood or saliva after cleaning needs to be prevented.
- Remove water residues with a faint airjet, but do not overdry dentine. The aim is to achieve a moist, but not wet dentine surface (Wet Bonding).

Colour design:

Dark or strongly yellowish necks of teeth impair the colour effect of the restoration. Therefore, these cervical areas should be masked with the enclosed **Opaker** system before applying the filling.

The set contains a shade guide for reference to show the colour effect of the **Opakers** combined with the restorative. The shade guide contains the following samples made from cured material:

- white - for a light gingival shade
- light - for a medium gingival shade
- dark - for a dark gingival shade
- nature - the basic gingival shade

The shade can be individually adapted step by step by mixing the **Opakers**.

Bonding material:

Amaris Gingiva is applied with a dentine enamel bond, using the adhesive technique. Any light-curing bonding materials can be used. Observe the relevant instructions for use during preparation (etching technique) and processing.

Light-curing:

Conventional polymerisation devices are suitable for light-curing the material. The light output should not be lower than 500 mW/cm².

Application of the Opaker:

Apply some **Opaker** onto the tooth surfaces to be masked and disperse thinly, then light-cure for **40 s**. If necessary, apply a second coat of **Opaker** to cover the tooth surface, and light-cure for **40 s**. In order to ensure the functionality of the non-dripping syringe do not pull back the plunger during or after use.

The attached application tip is intended for single use only. Remove it after use and seal the syringe again tightly with its original cap. The use of a protective barrier is recommended to protect the syringe against contamination with bodily fluids, dirt from hands or oral tissue.

Application of Amaris Gingiva:

Apply fillings of more than 2 mm thickness in layers. Light-cure each layer with blue light for **40 s**. Use transparent strips and light-wedges. Ensure a good marginal seal. Bring the light aperture as close as possible to the light-curing unit as otherwise the curing will be less effective. Inadequate curing can cause discolouration and complaints.

Finishing:

Contouring and finishing of the filling can be conducted immediately after removal of the shaping device, using cooling (e.g. fine or extrafine diamond polishers, polishing discs). Finally, the margin of the filling or preferably the entire tooth should be fluoridated.

Warnings, precautionary measures:

- Phenolic substances, especially products containing eugenol and thymol, lead to curing failures in restorative composites. Consequently, the use of zinc oxide/eugenol cements or other materials containing eugenol should be avoided in combination with restorative composites.
- Our information and/or advice do not relieve you of the obligation of checking that the products supplied by us are suitable for their intended purpose.

Constituents (in descending order according to content):

Amaris Gingiva: Barium aluminium borosilicate glass, prepolymer, UDMA, BisGMA, TEGDMA, methacrylated silica, initiators, stabilisers, pigments

Amaris Gingiva Opaker: TEGDMA, BisGMA, UDMA, HEDMA, fumed silica, initiators, stabilisers, pigments

Storage instructions and application method:

Storage at 4°C-23°C. Seal the syringes immediately after removing the material in order to avoid exposure to light and subsequent polymerisation. Do not use after the expiry date.

Disposal:

Dispose of the product in accordance with local regulations.

Reporting obligation:

Serious events such as death, temporary or permanent serious deterioration of a patient's, user's or other person's health condition and a serious risk to public health that arise or could have arisen in association with the use of **Amaris Gingiva** must be reported to VOCO GmbH and the responsible authority.

Note:

The Summary of Safety and Clinical Performance of **Amaris Gingiva** can be found in the European database on medical devices (EUDAMED – <https://ec.europa.eu/tools/eudamed>). Detailed information can also be found at www.voco.dental.

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VC 60 AA1970 E1 0322 99 © by VOCO