Opallis

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For professional use only

Carefully read the information contained in this document before using the product. Keep this document for later reference, at least until the product is completely consumed or until there is no further interaction between this product and the last patient in whom it will be used.

Product description

Opallis is a radiopaque, light-curing composite used to restore anterior and posterior teeth, with fluorescent and opalescent properties, high polishability and various shade options and translucency levels. The composite has average particle size of 0.5 microns, total inorganic filler by weight between 70% and 80%, and by volume between 57% and 58%.

Commercial presentation

Check the product presentations available in your country.

Refill: 1 syringe with 1 composite shade.

Syringe with 2 g: DA4, DC2, DBL1, DBL2, DBL3, D-BLEACH, A0.5, B0.5, EA4, EB3, EC2, EC3, EBL1, EBL2, EBL3, E-Bleach H, E-Bleach M, E-Bleach L, T-Blue, T-Neutral, Opaque Pearl (OP) and Opaque White (OW).

Syringe with 4 g: DA1, DA2, DA3, DA3.5, DB1, DB2, EA1, EA2, EA3, EA3.5, EB1, EB2, A1 Body, A2 Body and A3 Body.

Basic Kit: 6 syringes.

Shades: EA2, EA3, EA3.5, DA2, DA3 and T-Neutral.

Intro Kit: 4 syringes, 1 Condac 37% and 1 Ambar.

Shades: EA1, EA2, EA3 and DA2.

Clinical Kit: 13 syringes.

EA1, EA2, EA3, EB2, DA1, DA2, DA3, DB2, D-Bleach, T-Blue, T-Neutral, E-Bleach H and Opaque Pearl (OP).

Single-dose capsules: package containing 20 capsules with 0.2 g each.

Basic composition

The composite features radiopaque silane-treated glass particles (70-80%), methacrylic monomers (15-25%), silicon dioxide (2-7%), photo-initiator composition (< 1%), stabilizers (< 1%) and pigments (< 1%).

Intended Users

Qualified dental professionals who have theoretical and practical knowledge of the use of dental products.

Intended Patient Population

Patients who need to undergo a direct restorative procedure.

Clinical Benefit

Restoration of the form and function of compromised dental structure. In addition to promoting harmony and aesthetic balance of the smile.

Intended use

Restorative materials used for the purpose of treatment, attenuation or compensation of an injury or a deficiency, with the objective of efficiently restoring the form (anatomy), function and aesthetics of the parts of the compromised tooth structure.

Product indications

The composite is indicated for use on permanent and deciduous teeth:

- Direct restorations in anterior and posterior teeth (classes I, II, III, IV, V and VI);
- Direct veneers with composites;
- Cementation of dental fragments;
- Core building;
- Tooth veneering;
- Diastema closure or reduction;
- Modification of the tooth shape (e.g. conoid teeth);
- Porcelain or composite repairs.

Opallis shade system

Opallis uses the VITA Classical A1-D4® (VITA Zahnfabrik) shade guide as reference for its enamel and dentin shades. At the same time, it also offers shades with high translucency or opacity and shades for whitened teeth. In addition, the Opallis shade system includes Opallis Body resins, which have the ability to mimic both enamel and dentin, allowing for restorations with unique opacity.

- The two translucent shades (T-Blue and T-Neutral) allow for creating thin layers on top of the other shades, creating aesthetic restorations that perfectly reproduce the tooth structure.

The following is a brief description of the application of translucent shades to enamel:

- T-Blue: highly translucent composite, restoring the opalescent effect to the tooth.
- T-Neutral: translucent composite that reproduces the palatal side of the enamel and the incisal third. Its translucency is slightly less than the shades described above.
- The shades A0.5 and B0.5 can be used in pediatric dentistry to reproduce details of the anatomy of deciduous teeth.

- The extra-opaque shades, such as Opaque White (OW) and Opaque Pearl (OP), recreate details of the incisal halo and mask darkened dentin substrate.
- The composites for bleached tooth enamel (EBL1, EBL2, EBL3, E-Bleach H, E-Bleach M, E-Bleach L) and bleached dentin (DBL1, DBL2, DBL3, D-Bleach) were developed for restoring bleached teeth and deciduous teeth. Six colors are offered to reproduce the enamel and four colors to reproduce the dentin. The enamel color with lowest (lightest) saturation is EBL1, and the color with highest (darkest) saturation is E-Bleach L. In the dentin colors, the one with lowest (lightest) saturation is DBL1, and the color with highest (darkest) saturation is D-Bleach.
- Opallis resins must be cured at each increment with a blue light emitting curing light (400-500 nm) with 450 mW/cm² minimum power. Between increments, do not remove the oxygen-inhibited layer, as this will interfere with adhesion of the increments. Table 1 shows all available Opallis shades, the amount of product contained in each syringe and the recommended light-curing time.
- The light curing time described in the table must be followed to ensure proper depth of cure, degree of conversion, shade stability and mechanical properties of the composite. Overexposure is not possible.

Table 1: Available shades, amount and light curing time.

CATEGORY	SHADES	AMOUNT	LIGHT CURING TIME*
Dentin	DA1	4 g	40 s
	DA2	4 g	40 s
	DA3	4 g	40 s
	DA3.5	4 g	40 s
	DA4	2 g	40 s
	DB1	4 g	40 s
	DB2	4 g	40 s
	DC2	2 g	40 s
	DBL1	2 g	40 s
	DBL2	2 g	40 s
	DBL3	2 g	40 s
	D-Bleach	2 g	40 s
	A0.5	2 g	60 s
	B0.5	2 g	60 s

Enamel	EA1	4 g	20 s
	EA2	4 g	20 s
	EA3	4 g	20 s
	EA3.5	4 g	20 s
	EA4	2 g	40 s
	EB1	4 g	20 s
	EB2	4 g	20 s
	EB3	2 g	40 s
	EC2	2 g	40 s
	EC3	2 g	40 s
	EBL1	2 g	20 s
	EBL2	2 g	20 s
	EBL3	2 g	20 s
	E-Bleach H	2 g	20 s
Translucent Effect	T-Blue	2 g	20 s
	T-Neutral		
Opaque Effect	Opaque Pearl (OP)	2 g	60 s
	Opaque White (OW)		
Body	A1	4 g	40 s
	A2		
	A3		

*Considering light curing equipment with power above 450 mW/cm² and wavelength between 400 and 500 nm.

Precautions

- This product is intended to be used only in accordance with its usage instructions;
- Avoid the use of eugenol-containing materials, such as lining materials and temporary restoratives, as they may impair proper curing of the composite;
- Avoid contact of the uncured composite with the skin, mucosa and eyes as it can have an irritating effect due to the presence of methacrylates. The use of all personal protective equipment (PPE) is recommended when handling the product, such as gloves, masks and protective eyewear;
- When restoring large cavities, lining materials (such as calcium hydroxide) should be used on the pulpar wall;
- The composite is light sensitive, so be sure to close the syringe properly after use. The composite increment must be used right after it is removed from the tube;
- Wear protective eyewear when using curing lights;
- The product must not be swallowed or aspirated. If so, seek medical attention;
- The product is contraindicated if a dry operative field is not obtainable or if the restorative technique cannot be followed.

Contraindications

- Do not use the product if the patient is allergic to any of the substances present in the composition. In case of allergic reactions, discontinue use of the product and, if necessary, seek medical attention.

Side effects

This product contains methacrylate monomers,

which may cause allergic irritation in patients sensitive to these substances.

Directions for use

Carefully read the Warnings, Precautions, Contraindications and possible Side Effects in this document before using the product.

Shade selection

- Before selecting the shade, the tooth to be restored must be clean and hydrated, before isolation. Prophylaxis is recommended for this step;
- The VITA Classical A1-D4® (VITA Zahnfabrik) or OPALLIS® (FGM) shade guides should be used for shade selection;
- Alternatively, a small portion of the composite can be brought to the surface of the tooth to be restored (without prior etching or adhesive) and light cured for 30 seconds, then moistened with the patient's saliva and observed. This step can be performed to confirm the shade of the chosen composite.

Product application

1. Isolate the operative field. Isolation should ensure a moisture-free field with good visualization of the site. If necessary, use absolute isolation.
2. Prepare the cavity according to the adhesive technique.
3. Perform acid etching (enamel and dentin). Then wash with plenty of water, removing all acid from the surface, and dry the enamel gently, leaving the dentin moist.
4. Apply adhesive according to the manufacturer's recommendations.
5. Inserting the composite and light curing: the incremental technique is recommended because it

allows for better performance of the material. Increments of up to 1.5 mm must be inserted, shaped and light-cured (for the recommended time for each specific shade) in the cavity.

6. Fill the cavity in increments until the restoration is complete.

7. For finishing, use sandpaper discs or fine grit drill bits. For polishing, the use of felt discs and polishing pastes is recommended.

Note: Checking occlusion is extremely important. Overdenture restorations cause excessive stress on the tooth and the restoration, which can lead to pain. Do not leave restorations in infraocclusion.

Single-dose capsule:

1. Insert the capsule into the syringe/applicator gun system properly, so that it is ready for use.

2. Remove the cap from the capsule and apply the composite into the cavity by hand pressure on the syringe/applicator gun.

3. After completing the restoration (incrementally), discard the capsule.

Note: Opallis capsules are for single use only. Do not reuse, to avoid the risk of cross contamination.

Preservation and storage

- The product must be kept in its original packaging, properly closed;
- Protect the product from direct sunlight;
- Store the product at temperatures between 5 °C and 30 °C (41 °F and 86 °F);
- Do not freeze the product.

Note: If the product is stored in a refrigerator, preferably wait about 15 minutes before use. The best temperature range for applying the product is

20 °C – 25 °C.

Warnings

- This material was manufactured for dental use only and must be handled in accordance with the usage instructions;
- Do not use the product after its expiration;
- Follow the legislation in your country when disposing of the product;
- Keep out of the reach of children.

Expiration date

Printed on the product packaging.

This material has been manufactured for dental use only and must only be handled and stored in accordance with its usage instructions. The manufacturer is not liable for any damage caused as a result of use not in accordance with the guidelines, improper storage or incorrect handling. Moreover, the user is obliged to verify and confirm, before handling and/or applying the product, and under his/her responsibility, whether this material is compatible with the intended use. Furthermore, we underscore that data descriptions do not constitute any type of warranty and are therefore non-binding.

Legendas de símbolos Leyendas de los símbolos Explanation of symbols	
	Data de fabricação Fecha de fabricación Date of manufacture
	Número do lote Número del lote Batch Number
	Prazo de validade Fecha de validez Shelf life
	Limites de temperatura Límites de temperatura Temperature limits
	Não utilizar se a embalagem estiver danificada No lo utilice si el empaque está dañado Do not use if packaging is damaged
	Conservar seco Consérvelo seco Keep it dry
	Manter afastado da luz solar Manténgalo lejos de la luz solar Keep protected from sunlight
	Fabricante Fabricante Manufacturer
	Consulte as instruções de utilização Consulte las instrucciones de utilización Refer to instructions for use
	Representante na Comunidade Europeia Representante en la Comunidad Europea Representative in the European Community

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EN: Any serious incident that has occurred in connection with the device must be reported to the manufacturer and to the competent authority where the user/patient is based.
The summary of safety and clinical performance is available at: ifu.fgm.ind.br/sscp; just search by product name.