VOCO

Rebilda® DC

CE 0482

EN

Instructions for use

MD EU Medical device

In accordance with DIN EN ISO 4049

Product description:

Rebilda DC is a dual-curing, highly radiopaque flowable composite with excellent mechanical properties for core-build up.

Rebilda DC is available in three different shades to cover a large range of indications:

- blue to visualise the transition between material and tooth substance, especially in molar and pre-molar areas
- dentine shade for high aesthetic demands, e.g. under full ceramic restorations or composite restorations with high translucency
- white as a possibility to visualise preparation margins and simultaneously meet aesthetic demands

Indications:

- Adhesive core build-up of vital and non-vital teeth
- Luting of fibre-reinforced resin posts

Contraindications:

Rebilda DC contains methacrylates, BHT, amines and benzoyl peroxide.

Rebilda DC should therefore not be used in patients with a known hypersensitivity (allergy) to these ingredients.

Patient target group:

Rebilda DC is suitable for use in all patients without any age or gender restrictions.

Performance features:

The product's performance features satisfy the requirements of the intended use and the relevant product standards.

User:

Rebilda DC should only be used by a professionally trained dental practitioner.

Use:

- A) Take appropriate measures in areas close to the pulp.

B) Bonding material

Rebilda DC is used in the adhesive technique in combination with a light- or dual-curing adhesive which must be compatible with dual-curing composites. The product's instructions for use should be consulted to determine the correct preparation (etching technique) and procedure.

C1) Rebilda DC in the cartridge

Insert **Rebilda DC** cartridge into the **VOCO** dispenser (type 3) or any other extruder gun with a suitable piston. Due to technical reasons, there may be slight differences in the fill levels of both chambers prior to initial use. You should therefore remove the cap and continue to press out the material until an equal amount has been squeezed out of both outlets. Then insert a mixing tip type 12 and lock with a 90° turn clockwise. If necessary, attach an additional intraoral tip type 3 onto the mixing tip. For application in the root canal, use a mixing tip type 19 with integrated metal cannula. The cannula can be carefully pre-bent by hand. By pressing the lever the material is extruded and automatically mixed in the correct ratio.

C2) Rebilda DC in the QuickMix syringe

Due to technical reasons, there may be slight differences in the fill levels of both chambers prior to initial use. You should therefore remove the cap and continue to press out the material until an equal amount has been squeezed out of both outlets. Afterwards place a mixing tip type 11. The notches on the syringe and the mixing tip must line up exactly. Match the nooks on the syringe and the mixing tip and fix the mixing tip with a 90° turn clockwise. If necessary, put an application tip type 4 on top of the mixing tip. The material is mixed automatically in the correct ratio when dispensing and can be applied directly. All QM syringes come with a detachable piston. The piston can simply be pulled off from the QM syringe and reinserted. The plugs in the syringe chambers will stay in place preventing any cross contamination of the material in the syringe.

D) Application of the core build-up material

Place a matrix around the prepared tooth and apply **Rebilda DC** directly from the cartridge/QuickMix syringe.

Light-curing

The use of transparent matrices is advantageous. Place the light emission tip of the light-curing device (polymerisation light) as close as possible to the surface of the core build-up and polymerise for 40 s using an LED/ halogen lamp with a light output of at least 500 mW/cm².

If the distance is more than 5 mm, the curing depth may be reduced. In case of a larger core build-up, generally apply in layers of 2 mm and cure each layer for 40 s. After light-curing, the material can be finished immediately.

Chemical curing

Without light-curing, the build-up can be finished 5 min after its application in the mouth.

Rebilda DC has a working time of approx. 2 min.

In case a reshaped core form is used, apply **Rebilda DC** directly into the reshaped core form. Then seat the crown on the prepared tooth, adjust its position and light-cure the margins for 40 s to fix the piece.

The warranted product characteristics can only be achieved when using the supplied original mixing tips, which are also available as accessories. The attached mixing tip is intended for single use only. After use, store the cartridge/QM syringe tightly sealed. Check that material flows freely through the openings before reusing.

Luting of fibre-reinforced resin posts (in combination with a dual-curing adhesive):

1. When luting root posts, apply the dual-curing adhesive according to its instructions for use.

2. Using the cartridge:

When using mixing tip type 19, apply **Rebilda DC** directly into the root canal with the bendable metal cannula of the mixing tip. If necessary, the metal cannula of the mixing tip can be pre-bent by hand for this purpose. When not using mixing tips type 19, proceed as follows:

After preparing the post according to the manufacturer's instructions, wet post with **Rebilda DC**. Do not apply **Rebilda DC** directly into the root canal nor use a lentulo! Insert the post making sure to leave a slight excess.

Using QuickMix syringe:

It is recommended to apply **Rebilda DC** directly into the root canal using thin intraoral tips type 1. Do not use a lentulo! After preparing the root pin according to the manufacturer's instructions, it can be additionally moistened with **Rebilda DC**. Insert the post making sure to leave a slight excess.

3. To fix the post and cure the dual-curing adhesive applied to the occlusal surface, light-cure for 40 s.

4. For core build-up, follow the instructions given in D) Application of the core build-up material.

Warnings, precautionary measures:

- Phenolic substances, especially products containing eugenol and thymol, lead to curing failures in composites. Consequently, use of zinc oxide eugenol cements or other materials containing eugenol must be avoided.
- Do not use bonding agents which are not compatible with chemically curing or dual-curing composites.
- Insulate the core build-up before taking impressions or making temporary crowns/bridges (good moistening is often sufficient) to avoid bonding to these materials.
- Avoid contact with the oral mucosa.
- Our information and/or advice do not relieve you of the obligation of checking that the products supplied by us are suitable for their intended purpose.

Constituents (in descending order according to content):

Barium aluminium borosilicate glass, UDMA, HEDMA, DDDMA, pyrogenic silica, BisGMA, initiators, stabilisers, pigments

Storage:

Do not expose **Rebilda DC** to direct sunlight or the chairlight. Store at 4°C–23°C. Do not use after the expiry date.

Disposal:

Dispose of the product in accordance with local regulations.

Reporting obligation:

Serious events such as death, temporary or permanent serious deterioration of a patient's, user's or other person's health condition and a serious risk to public health that arise or could have arisen in association with the use of **Rebilda DC** must be reported to VOCO GmbH and the responsible authority.

Note:

The Summary of Safety and Clinical Performance of **Rebilda DC** can be found in the European database on medical devices (EUDAMED – <https://ec.europa.eu/tools/eudamed>).

Detailed information can also be found at www.voco.dental.

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VC 60 BB1395 E1 0622 99 © by VOCO