Whiteness HP Blue

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INSTRUCTIONS FOR USE

PT	Português Clareador Dental para Dentes Vitais	02
ES	Español Aclarador Dental para Dientes Vitales	06
EN	English Whitening gel for Vital Teeth	10

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Whitening gel for Vital Teeth

For Professional Use Only

Read the information in this manual carefully before using the product. Save it for reference at least until the product is completely used and/or until there is no more interaction between the product and the last patient.

Product Description

Whiteness HP Blue is a whitening gel with a base of 35% hydrogen peroxide in its formula for dental in-office use. It is sold in two phases contained in syringes that must be attached and mixed for the application. The application of the gel involves a single procedure during the session and a violet pigment which makes visual control easier and the color does not change during the process. Whiteness HP Blue has calcium in its formula and its pH is stable and alkaline, minimizing enamel demineralization as a result of the whitening. Besides, the activating composition of the product makes it more active and efficient, dispensing the use of external sources of acceleration.

Product Indication

The product is recommended for the following cases:
Naturally dark or yellowish vital teeth.

Vital teeth darkened by age.

Vital teeth that are darkened or stained by extrinsic substances such as pigments from coffee, tobacco, tea, etc.

Vital teeth with light to moderate stains due to fluorosis.

Vital teeth stained by tetracycline, in case of light stains.

In more challenging cases (severe darkening) the

product can be combined with at-home whitening technique.

Basic composition

Active Ingredients: hydrogen peroxide at 35% (after mixing the two phases).

Inactive Ingredients: thickeners, inert pigment, neutralizing agents, calcium gluconate, glycol and purified water.

Commercial Availability

Check the commercial availability your country.

Whiteness HP Blue - Kit

06 sets of syringes (1,2g per kit, 7,2g total of gel).

01 bottle with 2g of Neutralize.

01 syringe of Top dam with 2g.

06 syringe attaching devices.

12 applicator tips.

Refill

01 set of syringes (1,2g total of gel).

01 syringe attaching device.

01 syringe of Top dam with 1g.

03 applicator tips.

Precautions

The product is for exclusive use in the office and must be handled by dental professionals.

The patient must be in good buccal and general health condition before beginning the treatment. During product handling, both the professional and the assistant must wear protective gloves and goggles. The patient must also wear protective goggles and other necessary protections to avoid accidental contact of the product with the skin and clothes.

The whitener is composed by a phase called peroxide

(phase 1) and another phase called thickener (phase 2), which must be mixed at the moment of use. The professionals must be careful when manipulating the peroxide phase because it contains hydrogen peroxide at 50%. The content of the syringes with peroxide may be under pressure and/or may have leaked due to adverse storage conditions. Before opening the package, protect hands and eyes (professional dentist and assistant) and observe if there is no liquid leaking in its interior. If there is a leak, request replacement of the material.

Connect the syringes by means of the attaching device that comes with the product and make sure they are securely attached. Failure to do so may lead to leaking of the product.

After mixing the phases, insert the tip and check to make sure there is no resistance for the product to flow before applying to the patient's teeth. If there is any resistance, replace the tip and test on another surface (example: a mixing pad).

Whiteness HP Blue is a whitener that is highly oxidizing and when in contact with live tissues may cause a white stain and temporary irritation (usually disappearing in 2 hours or less, without further complications).

In case of accidental contact with the skin or intraoral live tissues, inhibit the peroxide covering the area with a few drops of Neutralize solution and wash.

Patients must be oriented to avoid ingesting acidic and/or highly colored foods for at least 24 hours after the whitening process to avoid sensitivity and not to occasionally compromise the treatment.

Inform the patient about the Pros and Cons of the whitening procedure.

The product must be used according to the recommendations in the Instructions for Use.

Contraindications

It is not necessary to use light sources (LED, high-intensity LASERs or other systems), since those light sources do not improve results and do not reduce the gel application time.

Dental enamel acid etching previous to whitening is not necessary and neither recommended.

The product is not recommended for patients who have recently undergone gingival surgery or that have periodontal diseases, sensitive teeth, exposed roots, cavities, ulcers or bad health in general. In those cases, the solution of the problem must be provided before beginning the treatment.

The product is not recommended for whitening teeth with amylogenesis and dentinogenesis imperfection, severe fluorosis, intense tetracycline stains and other anomalies of the enamel and dentin that may compromise the vitality of the teeth.

The product is not recommended for the use in patients under anesthesia. It is important to monitor the sensitivity felt by the patient during the whitening. In extreme cases in which there is an imperfection or flaw in the dental structure that allows for the peroxide to reach the dental pulp in greater amounts, pulp death may occur.

The product is not recommended for use in pregnant or lactating women or in patients younger than 15.

The product is not recommended for patients allergic to the components of the product.

Side effects

Dental hypersensitivity may occur during or after whitening. Since this is hard to predict in several situations, the use of the desensitizer Desensibilize KF 0.2% - at -home use - (FGM) or Desensibilize KF

2% (FGM) – in-office use - is recommended for a prophylactic desensitization before the beginning of the treatment.

Dental whitening may increase the visibility of pre-existing stains due to hypoplasia of the enamel (opaque white stains). In most cases, no additional care is necessary. The stains will disappear after the end of the treatment and the teeth will recover their natural hydration. For more resistant cases, treatment with fluoride in the area may help its disappearance.

Mode of use

Preparing for Whitening

- Put on (dentist and assistant) your safety equipment (gloves, goggles, cap, mask and scrubs) to handle the product. Also, make sure your patient wears protection to avoid accidental contact with the gel and an adequate lip retractor (Arcflex – FGM).

- Remove the whitener from the package following the safety recommendations.

- Remove the cap from the syringe with the thickener and attach the device between the syringes.

- Carefully attach the syringe with the peroxide noting that its content may be under pressure.

- Mix the content of the two syringes only at the moment of application.

Step by step

Whitening Vital Teeth



Fig 1. Position the ArcFlex (FGM) lip retractor for retracting lips and cheeks.

Fig 2. Carry out the prophylaxis with pumice stone and water.

Fig 3. Register the initial shade of the teeth with a shade guide and/or photograph the initial case.

Fig 4. Apply Desensibilize KF 2% (FGM) for 10 minutes to promote a previous desensitization. Remove the excess gel with a suction tip, wash and dry.

Fig 5. Promote relative isolation with the light cured gingival barrier Top dam (FGM) covering the marginal gum and the papillae with a layer of 3 to 5mm of width and 1-1.5mm of thickness. The protection must cover approximately 0.5 to 1mm of the dental surface. Check if the product is covering the whole surface of the gums.

Fig 6. Light cure the Top dam composite (FGM) for 20 to 30 seconds for each group of 3 teeth. The gingival protector that forms is rigid and insoluble, preventing possible irritation by aggressive products.

Fig 7. Mix the two phases with the connected syringes, pushing the drums alternately for at least 8 times (4 times on each side), and then push the whole mixed content to one of the syringes. Then it is ready for use.

Fig 8. Attach one tip to the syringe that holds the gel and apply one layer over the whole vestibular face of the teeth to be whitened (including the interproximal)

and extend a little towards the incisal and occlusal faces. The gel layer should have between 0.5 and 1mm of thickness. Consider the smile line (generally from the second premolar to the second premolar of the other arch).

Note: a kit of syringes with 1.2g of gel is enough for one simultaneous application on both arches in the smile line.

Fig 9. Leave the gel on the dental surface for up to 40 minutes. With the help of a Cavigbrush micro-applicator move the gel over the teeth frequently (every 5 to 10 minutes) to release possible oxygen bubbles and renew the contact between the gel and the teeth.

Note: the duration of application of whiteners are suggestions and should be managed by the professional. Shorter times or interruptions of the treatment may be used in cases of dental hypersensitivity.

Fig 10. At the end of the treatment, remove the gel with an endodontic or surgical cannula and wash the teeth with abundant water. Remove the gingival barrier by detaching it with an explorer probe.

Fig 11. Polish the teeth with Diamond Excel polishing paste and Diamond or Diamond Flex felt.

Important notes

1. Monitor the patient during the entire duration of the whitening treatment for dental sensitivity and possible irritated sites caused by contact with the peroxide. In case any of this occurs, suspend the process and analyze the possibility of an intervention (refer to the comments in the section Warning and Side Effects).

2. Most of dental stains may be removed with only one whitening session. However, depending on the origin and the age of the staining, there are cases in which more than one session may be necessary.

3. If the whitening obtained in one session does not fulfill the expectations and the patient does not experience sensitivity or any other contraindication, the reapplication of the product may be carried out for up to two more sessions. In case more than one session is necessary, a minimum interval of 7 days between sessions should be observed.

4. In case the patient develops sensitivity that cannot be controlled, the professional may alternatively decrease the whitening times and divide it between 3 sessions. In that case, it is suggested that the professional check if there is no abnormality in the dental structure that may be causing the sensitivity (fissures, dentin exposure, etc.). When the sensitivity is difficult to control or, for any reason, a new whitening session is not possible or recommended, the whitening treatment may be complemented with at home use whiteners (Whiteness Perfect at 10% of carbamide peroxide, for example).

5. Like any other whitening technique, the shade record of the teeth before treatment, the warnings to the patient about the real perspectives of his/her case (limitations in cases of tetracycline stains, greyish color, etc.), warnings about possible sensitivity, need to replace restorations, are recommended. There are cases in which, due to particular characteristics of the teeth (enamel, type of pigmentation, etc.) the whitening level desired may not be possible.

Conservation and Storage

Store the product under temperatures between 5°C and 25°C / 41°F and 77°F.

Do not freeze the product.

Keep away from direct sunlight.

Handling conditions

During product handling, both the dentist and the

assistant must wear protective gloves and goggles. The patient must also wear protective goggles and other necessary protections to avoid accidental contact of the product with the skin or clothes.

Handle the product away from the patient over a sink or in a safe space and without risks of contaminating other people.

After mixing the whitener, it must be used in the same session. The gel may not be stored for use for a later session after it is handled, as it will lose its efficacy.

In case there is any gel remaining, it must be disposed (diluted in water and disposed in a sink, for example).

For product disposal, check your country's disposal legislation.

Warnings

Do not use the product after its expiration date.

Keep out of the reach of children.

This material is for dental purposes only and should be handled and stored only in accordance with its instructions for use. The manufacturer will not be held responsible for any damage caused as a result of misuse, improper storage or improper handling. In addition, the user is required to check and verify, before handling and / or application of the product, and under their responsibility, whether this material is compatible with its intended use. Moreover, it is emphasized that data descriptions do not constitute any kind of guarantee and, therefore, have no binding.

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