
# Unimatrix® Hard ○ Soft

- PT Sistema de matrizes seccionais pré-formadas para restaurações de Classe II
- ES Sistema de matrices seccionales pre-formadas para restauraciones Classe II
- EN Preformed sectional matrices system for Class II restorations
- FR Système de matrices sectionnelles pour des restaurations de Classe II
- DE System für vorgeformte Teilmatrizen für Restaurierungen der Klasse II
- IT Sistema di matrice sezionale preformata per i restauri di Classe II
- NL Systeem van voorgevormde sectionele mallen voor klasse II restauraties
- CS Systém práce s sekčními matricemi u oprav třídy II

Hard

Soft

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ES

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**RING / SILICONE TIPS - Cleaning:** Carry out the cleaning with the aid of water and enzymatic detergent 4 enzymes (Amylase, Protease, Lipase and Carbohydrazase) and with pure pH range between 6.5 - 7.5. Prepare the solution according to the manufacturer's instructions. The water should be potable. For a better effectiveness, it is recommended that the temperature be between 40°C and 55°C. Submerge the products and leave in contact with the solution. Remove parts and rinse under running water for at least 10 seconds. Dry the pieces with absorbent paper or disposable wipes. **Sterilization:** The ring and the triangular silicone tips must be sterilized separately in autoclave until 134 °C - 6 min. The sterilization can cause slight changes in the coloring of the tips and reinforcement, but does not affect its functionality. Check the product before each use. Note if there is damage (breaks/cracks/rust) that prevents proper operation and does not use in such cases. These conditions are indicative of fragility of the material, whose strength can be reduced during clinical use. **Retreatment Limits:** This is a reusable device. Repeated processing has minimal effect due to its composition. The end of the life of the ring and triangular silicone tips, identified by the professional, is determined due to damage resulting from use.

**Composition:** Matrices: stainless steel. Ring: nickel-plated carbon steel. Ring reinforcement: thermoplastic. Tips: silicone

**Storage:** To be stored in a cool and dry place, protected from sunlight.

**Disposal:** Dispose in accordance with applicable local requirements.

**Ref.** 004126H  
004124R  
004122R  
004120R  
004126S  
003020102  
004159T  
004160T  
-  
**Code** 0003020081  
0003020085  
0003020086  
0003020087  
0003020091  
0003020103  
0003020104  
0003020105  
0003020106

Unimatrix R anillo Hard – 1 anillo + 4 protectores de silicona triangulares  
Unimatrix R – 50 matrices (P/M/G) + 2 anillos + 8 protectores de silicona triangulares  
Unimatrix R – 25 matrices (P/M/G) + 1 anillo + 4 protectores de silicona triangulares  
Unimatrix R – 100 matrices (P/M/G) + 2 anillos + 8 protectores de silicona triangulares  
Unimatrix R ring Soft – 1 ring + 4 triangular silicone tips

**Ref.** 004126H  
004124R  
004122R  
004120R  
004126S  
003020102  
004159T  
004160T  
-  
**Code** 0003020081  
0003020085  
0003020086  
0003020087  
0003020088  
0003020103  
0003020104  
0003020105  
0003020106

Unimatrix R – 50 matrices (S/M/L) + 2 rings + 8 triangular silicone tips  
Unimatrix R – 25 matrices (S/M/L) + 1 ring + 4 triangular silicone tips  
Unimatrix R – 100 matrices (S/M/L) + 2 rings + 8 triangular silicone tips

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**Instructions for use** EN

The Unimatrix R System is indicated for Class II restorations and composed by preformed sectional matrices, available in 3 sizes, and reinforced rings. The main advantages of this system are the easiness in the reconstruction of contact points and natural tooth anatomy. The reinforced rings are equipped with triangular silicone tips, and available in two models: Hard (blue) indicated for clinical cases requiring more pressure and dental separation, and Soft (pink) indicated for clinical cases requiring less pressure and dental separation.

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